

## Liquid Argon Technology for Neutrinos

Mitch Soderberg 7th IDS-NF Plenary Meeting





#### Talk Outline

• Introduction to Liquid Argon detectors for neutrinos.

The MicroBooNE Experiment

Future directions.

#### Introduction

- Liquid Argon Time Projection Chambers (LArTPCs) combine finegrained tracking and calorimetry, and appear scalable to very large size.
- •U.S. efforts to develop LArTPCs have expanded significantly in recent years.
- These efforts are aimed at developing the technology for a multi-kiloton detector that could be used to do a variety of physics (accelerator neutrinos, proton decay, astrophysics, ...)
- "Smaller" scale detectors, O(100 tons), also offer opportunity for important physics measurements.


## Why Noble Liquids for Neutrinos?

- Abundant ionization electrons and scintillation light can both be used for detection.
- If liquids are highly purified (<0.1ppb), ionization can be drifted over long distances.
- Excellent dielectric properties accommodate very large voltages.
- Noble liquids are dense, so they make a good target for neutrinos.
- Argon is relatively cheap and easy to obtain (1% of atmosphere).
- Drawbacks?...no free protons...nuclear effects unavoidable.

|                             | 9      | Ne     | Ar     | Kr     | Xe     | Water |
|-----------------------------|--------|--------|--------|--------|--------|-------|
| Boiling Point [K] @<br>1atm | 4.2    | 27.1   | 87.3   | 120.0  | 165.0  | 373   |
| Density [g/cm³]             | 0.125  | 1.2    | 1.4    | 2.4    | 3.0    | 1     |
| Radiation Length [cm]       | 755.2  | 24.0   | 14.0   | 4.9    | 2.8    | 36.1  |
| dE/dx [MeV/cm]              | 0.24   | 1.4    | 2.1    | 3.0    | 3.8    | 1.9   |
| Scintillation [γ/MeV]       | 19,000 | 30,000 | 40,000 | 25,000 | 42,000 |       |
| Scintillation λ [nm]        | 80     | 78     | 128    | 150    | 175    |       |

## Liquid Argon Neutrino Detectors


- Neutrino interactions in the TPC produce charged particles that ionize the argon as they travel.
- Ionization is drifted along E-field to wireplanes, consisting of wires spaced ~millimeters apart.
- Location of wires within a plane provides position measurements...multiple planes give independent views.
- Timing of wire pulse information is combined with known drift speed to determine drift-direction coordinate.



Images from ICARUS 50-liter TPC.

## Liquid Argon Efforts at Fermilab

Development focused on scaling LArTPCs to sizes necessary for long-baseline experiment.




#### Refs:

<sup>1.)</sup> A Regnerable Filter for Liquid Argon Purification Curioni et al, NIM A605:306-311 (2009)


<sup>2.)</sup> A system to test the effect of materials on electron drift lifetime in liquid argon and the effect of water Andrews et al, NIM A608:251-258 (2009)

## Liquid Argon Worldwide

#### ICARUS Program




#### Europe/Japan Program




direct proof of long drift path up to 5 m








Charged particles test beam, calorimetry, non-evacuated vessels, LAr purity Full
engineering
demonstrator
for larger
detectors, with
a stand-alone
short baseline
physics
programme




31










### MicroBooNE

- MicroBooNE will operate in the Booster neutrino beam at Fermilab starting in late 2013.
- Combines timely physics with hardware R&D necessary for the evolution of LArTPCs.
  - MiniBooNE low-energy excess
  - ▶ Low-Energy Cross-Sections
  - ▶ Cold Electronics (preamps in liquid)
  - Long drift (2.5m)

| Cryostat Volume               | 150 Tons                         |  |
|-------------------------------|----------------------------------|--|
| TPC Volume (l x w x h)        | 89 Tons (10.4m x 2.5m x 2.3m)    |  |
| # Electronic Channels         | 8256                             |  |
| Electronics Style (Temp.)     | CMOS (87 K)                      |  |
| Wire Pitch (Plane Separation) | 3 mm (3mm)                       |  |
| Max. Drift Length (Time)      | 2.5m (1.5ms)                     |  |
| Wire Properties               | 0.15mm diameter SS, Cu/Au plated |  |
| Light Collection              | ~30 8" Hamamatsu PMTs            |  |

- ★Stage 1 approval from Fermilab directorate in June 2008
- ★DOE CD-0 (Mission Need) in October 2009
- **★**DOE CD-1 June 2010
- ★DOE CD-2/3a (September 2011)



- → Joint NSF/DOE Project
- →\$1.1M NSF MRI for TPC, PMTs

#### MicroBooNE: Location

- MicroBooNE will sit on surface in on-axis Booster beam (BNB), and off-axis NuMI beam.
- Liquid Argon Test Facility will be located directly upstream of MiniBooNE enclosure. L=470m.
- Large event samples will allow a variety of cross-section measurements.

|                    | BNB         | NuMI        |
|--------------------|-------------|-------------|
| Total Events       | 145k        | 60k         |
| $\nu_{\mu}CCQE$    | 68k         | 25k         |
| NC πº              | 8k          | 3k          |
| $\nu_{\rm e}$ CCQE | 0.4k        | 1.2k        |
| РОТ                | $6x10^{20}$ | $8x10^{20}$ |


Projected Event Rates for MicroBooNE in 2-3 years.



Neutrino Beams at Fermilab



MicroBooNE will be located in Liquid Argon Test Facility




- L010/185kA
- ME

10<sup>-3</sup>
10<sup>-2</sup>
10<sup>-3</sup>
0.0 0.5 1.0 1.5 2.0 2.5 3.0 3.5 4.0 4.5 5.0 ν<sub>μ</sub> Energy (GeV)
NuMI Off-Axis Flux

## MicroBooNE: Physics

- Address the MiniBooNE low energy excess
  - MiniBoone is a Cerenkov detector that looked for  $v_e$  appearance from a beam of  $v_{\mu}$
  - Does MicroBooNE confirm the excess?
  - Is the excess due to a electron-like or gamma-like process?

#### MiniBooNE $v_e$ Appearance Result



MiniBooNE Result Excess

200-300MeV: 45.2±26.0 events

300-475MeV: 83.7±24.5 events


MicroBooNE will have ~5.5σ significance for electron-like excess, ~4σ for photon-like excess.

#### Refs:


1.) Unexplained Excess of Electron-Like Events From a 1-GeV Neutrino Beam MiniBooNE Collaboration, Phys. Rev. Lett. 102, 101802 (2009)

## MicroBooNE: Physics

- Prove effectiveness of electron/gamma separation technique (using dE/dX information).
- Low Energy Cross-Section Measurements (CCQE, NC  $\pi^{o}$ ,  $\Delta \rightarrow N\gamma$ , Photonuclear, ...)
- Continue development of automated reconstruction (building on ArgoNeuT's effort).




dE/dx in first 5cm of Simulated/Reconstructed Gamma/Electron showers.




Inclusive CC cross-section in neutrino-mode (simulated/ reconstructed/analyzed MC)

#### MicroBooNE: TPC

- TPC has 3 instrumented wireplanes (Two Induction at +/-60 from vertical, One Collection with vertical wires).
- Cathode is held at -125kV, setting up 500V/cm drift field.
- Wires are individually terminated around brass ferrules, then positioned on wire carriers.




Schematic of MicroBooNE TPC



Prototype wires and wire carrier boards.


#### MicroBooNE: Electronics

- CMOS preamplifiers located in liquid, attached to TPC.
- 12-bit ADCs sampled at 2MHz (i.e. 500ns per sample) for 4.8ms (x3 drift window).
- 1-hour data buffering for Supernova detection signal from SNEWS.



## MicroBooNE: Cryogenics

- Cryogenic system consists of filters/pumps/etc... for circulating and purifying LAr.
- Cryostat is evacuable (though the plan is not to evacuate) and foam insulated.




Schematic of MicroBooNE Layout




LAPD @ Fermilab

### Massive LArTPC Detectors

- Description here is the Reference design for the LBNE project.
- LArTPC at DUSEL would be two ~20 kTon modules.
- Detector located at the 800ft level at DUSEL.



~20 kTon LArTPC module(s)



800-ft. level layout.

| Cryostat Volume           | ~25 kTons           |  |
|---------------------------|---------------------|--|
| TPC Volume                | ~16.7 kTons         |  |
| # Readout Wires           | ~645000 (128:1 MUX) |  |
| Wire Pitch                | ~3 mm               |  |
| Electronics Style (Temp.) | CMOS (87 K)         |  |
| Max. Drift Length         | ~2.5m               |  |
| Light Collection          | TBD                 |  |

# Liquid Argon for Neutrino Factory

## Conclusion

- Liquid Argon detectors provide exceptional capabilities for neutrino physics, and there is significant R&D ongoing at Fermilab, and worldwide, to develop this technique for very large scales.
- MicroBooNE is the next major step in the U.S. plan for LArTPCs, and it will provide interesting physics and hardware development.
- Neutrino Factory...