






# A New Event Display Analysis Tool for MicroBooNE Handscan

---

Ellen Klein  
Yale University

MicroBooNE Summer Student  
Presentations  
08/03/12


# Introduction to TPC




ArgoNeuT Data!

# LArSoft Reconstruction

## Reconstruction Chain for LArTPC Data



Wire Readout



Reconstructed Hits

E. Klein





# Goal

---

- Create hand scan tool where evd users can select hits in a cluster (add/remove individual hits), run reconstruction and calorimetry algorithms on clusters, look at calorimetry/PID
- Evd will create clusters, all other reconstruction algorithms should work as is
  - Change fcl parameter so that use evd clusters


# Hit Selection

- Draw lines over particle tracks in evd
  - Direction is important!  
Keeps track of start and end points
- Evd will select hits within a 1.5 cm box around line, selected hits are blue
- Box width is an option set by user





# Hit Selection (2)

- Shift+click to add or remove hits from the cluster



# Reconstruction

- Clicking “Reload” button will pass selected hits to clustering producer




- Can pass these clusters to other reconstruction algorithms: track3Dreco, spacepts, spacepoints, seeds, etc.

Thanks to Tingjun and Ornella  
for calorimetry and PID code  
and Mitch for display


# Reconstruction (2)

Can also run calorimetry/PID on clusters from the  
evd to help determine particle type!






# Muon Reconstruction




# Data



# MicroBooNE

- Works in MicroBooNE geometry too!
- Note clear resolution on delta ray
- Still exploring which tracking algorithms are best
- Plan to use this for MicroBooNE oscillation handscan





# Current Status

---


- Fixing one last issue, then will check in the new Event Display code
- Evd creates clusters, all other reconstruction algorithms work without changes
  - Change fcl parameter so that they use evd clusters
- Step by step instructions for use are in backup slides



## Next Steps

---

- Test with other reconstruction algorithms (ie seeds)
- Create multiple tracks per event and save reconstruction to event record (shouldn't matter for hand scan)
- Any other ideas? Let me know!  
(Ellen.Klein@Yale.edu)



Thank You!




# Backup

---

# Extras

- Box on left side of evd will display track parameters (length, phi, theta)
- Information calculated from the two lines that user draws NOT from track reconstruction







# Step by Step (1)

- Decide which reconstruction algorithms you want to run, add these to your evd fcl file as producers
- Make sure they will look for the clusters produced by the evd (graph)

```
physics:
{
  producers:
  {
    graph: @local::argoneut_graphcluster
    spacepts: @local::argoneut_spacepts
    calo: @local::argoneut_calorimeter
    pid: @local::argoneut_chi2pid
  }

  filters:{}


  analyzers:
  {
    evdisp:{module_type: EVD}
    angles: @local::standard_spacepointana
  }
}

#list the modules for this path, order matters, filters reject all following items
evd: [ evdisp ].
#recopath: [showeranglecluster, showerfinder]
recopath: [graph,spacepts, calo ,pid]
ana: [angles]
#end_path are things that do not modify art::Event, includes analyzers.
#and output modules. all items here can be run simultaneously
trigger_paths: [recopath]
end_paths: [ana,evd],.
}

physics.producers.showeranglecluster.ClusterModuleLabel: "graph"
physics.producers.spacepts.ClusterModuleLabel: "graph"
physics.producers.calo.TrackModuleLabel: "spacepts"
physics.producers.pid.TrackModuleLabel: "spacepts"
```


## Step by Step (2)

- Run the event display on a file that has reconstructed hits
- Zoom views so that you can clearly see the track you wish to reconstruct
  - “Zoom Interest” is helpful
- Check “Make Clusters”




## Step by Step (3)

- Click and drag mouse to draw a line over the track you want to reconstruct
  - Click on track vertex, then drag to end of track
- Click “Select Hits” to select hits in a box around the line, selected hits will turn blue
  - Adjust box width
  - Don’t have to redraw line, just click “Select Hits” again





# Step by Step (4)

- Shift+click to add/remove individual hits from cluster
  - Useful for overlapping tracks, delta rays, etc.
  - Helpful to zoom in on parts of the track, especially end points!
- “Clear Hits” button will unselect hits in ALL planes




# Step by Step (5)

- When satisfied with selected hits, click “Reload” to make clusters and run reconstruction modules



# Step by Step (6)

- Mitch added calorimetry window to the event display
- After reloading, click Window->Calorimetry
- Useful to look at particle ID to determine particle/event type



# Step by Step (7)

- Click “Next” to clear everything and go to next event in the run
  - Close the calorimetry window first, else will crash!  
-- will fix this though

