

The Laser calibration system for the MicroBooNE detector

A. Ereditato, I. Kreslo, M. Lüthi, C. Rudolf von Rohr, T. Strauss, M. Weber

Albert Einstein Center for Fundamental Physics,
Laboratory for High Energy Physics,
University of Bern

Outline

- > The “laser”
- > Principle of multiphoton ionization
- > Results of laser measurements
- > Feedthrough design
- > Feedthrough test facility
- > Beam line-up
- > Power supply, racks, interlocks
- > Next steps
- > Conclusion

The “laser”

- > Continuum Surlite I-10
- > Nd:YAG laser with 1064 nm output
- > Wavelength is reduced to 266 nm by two second-harmonics generators
- > Laser properties:
 - Max pulse frequency 10Hz
 - Energy max. 60 mJ
 - Pulse width 4-6 ns
 - Divergence 0.5 mrad



One laser was ordered last week

Multiphoton ionization


- > $\lambda = 266 \text{ nm}$ correspond to $E_\gamma = 4.67 \text{ eV}$
- > For ionization $\sim 14 \text{ eV}$ are needed
- > For non-resonant states the lifetime is given by

$$\tau_\gamma = \frac{\lambda}{2\pi c} = 1.4 \times 10^{-16} \text{ s}$$

- > For quasi-resonant states one has

$$\tau_\gamma \propto \frac{1}{\Delta E} = \frac{1}{(E_i - E_\gamma)}$$


- > The laser has to have enough intensity to allow a three-photon ionization



See also: B. Rossi et al. 2009 JINST **4** P07011

I. Badhrees et al. 2010 New J. Phys. **12** 113024



Results of laser measurements



This picture shows a laser event in the Argontube detector. In the upper display the collection view and below the induction view are shown.




Results of laser measurements

- > Recent measurements with Argontube prove the feasibility of field corrections
- > Above: collection view of a laser path in distorted field
- > Below: same event after the electric field correction is applied




Results of laser measurements

- > Right side: correction applied to a cosmic muon event
- > Bottom: Reconstruction of the Electric field in Argontube




Feedthrough design

- > Quartz glass tube immersed in liquid Argon, providing a flat surface for the laser
- > CF-flange with a turnable head in order to swivel the mirror horizontally
- > A linear feedthrough allows to tilt the mirror vertically
- > The flange has to be vacuum tight



Feedthrough design



Feedthrough design

Top mirror holder

Turnable part

Fix part

Linear feedthrough bellow



Glass tube holder

Bearing ring




Feedthrough test facility


- > Test will be performed in medium-Argontube
- > Cylindrical vessel
 - Height 1.1 m
 - Diameter 71.4 cm
- > Vacuum test of the feed-through down to 10^{-5} mbar
- > Torque and force measurements of the movable parts
- > Cold test in liquid Argon



Beam line-up




Beam line-up



View from the point where the cold mirror will be placed

Power supply, Racks, Interlocks

- > Input 220 V single phase for the laser
- > Interlock with PMT power supply (Laser off, if PMTs on)
- > Interlock of the power supply if laser box is opened
- > Panic button switches of laser if pushed
- > System is designed with help of Linda, Bryce and Matt Quinn (LSO)




Next steps

- > Assembly of the feedthrough prototype (Nov 2012)
- > Test of the feedthrough prototype (Nov 2012)
- > Laser interlock system design (Nov 2012)
- > Steering motor selection and order (Dec 2012)
- > Arrival of the first laser in Bern (Dec 2012)
- > Alignment test of the laser in Bern (Jan 2013)
- > Alignment tests at Fermilab (Apr 2013)

Conclusion

- > Argon can be ionized by an UV-laser
- > The laser system is able to map the electric field
- > Design phase of the feedthrough and the beam line-up is completed
- > Feedthrough prototype parts are manufactured and ready for assembly
- > Optical instruments are ordered
- > First laser will arrive in Bern in about 7 weeks
- > Preparations of feedthrough tests has started

Tank you!



Backup


u^b

b
**UNIVERSITÄT
BERN**

AEC
ALBERT EINSTEIN CENTER
FOR FUNDAMENTAL PHYSICS

Why a laser?

- > Cosmic rays produce Ar^+ -ions
- > Ion drift velocity is only of the order of cm/s
- > Ar^+ accumulates
- > Field distortions occur
- > Laser delivers a straight path
- > Field distortions can be corrected



Eric Voirin: MicroBooNE-doc-1895-v4

Laser position

